Can be controlled by all known Cardio Pulmonary Exercise devices







Highlights

Reliable and reproducible stress tests

The experience of professionals who calibrate many ergometers showed that the Lode ergometers are the most reliable across the complete workload and rpm range and still within specifications even after many years of intensive use.

High standards

Lode is a socially and environmentally responsible company. All Lode products are RoHS/WEE compliant and Lode is ISO 9001:2003, ISO 13485:2008 and FDA 510K certified. All medical products comply to MDD 93/42/EEC, incl. IEC 60601-1.

Various test modes

Besides the hyperbolic (rpm-independent) mode that is used most of the time, the standard control unit offers several other test modes, like the fixed torque mode and the linear mode. These modes can be used in both manual and terminal mode.

Q-factor equal to road-bike

The Q-factor of the ergometer is equal to the Q-factor of road bikes, creating perfect training circumstances.

Rotatable handlebar with new lever

The new designed lever makes it even easier to adjust the handlebar. The handlebar can be rotated 360 degrees and is constructed in such a way that the test subject can be installed comfortably at every seating height.





Can be controlled by all known Cardio Pulmonary Exercise devices



The Corival is one of the most popular ergometers worldwide. The low start-up load of 7 Watt is first-class. The Corival cpet is standard supplied with a communication module and can therefor be easily controlled by all known stress ECG and pulmonary devices in the world. The workload, rpm and time can be readout from the 3,5" colour display. The Corival has an eddy current electro-magnetic braking mechanism. The biggest advantage of this system is the accuracy which is one of the most important Lode principles. With this ergometer, the stress tests performed are reliable and reproducible. The workload is adjustable in a range of 7 to 1000 watt. The ultralow step-through enables easy access to the ergometer and the latest design guarantees a perfect ergonomic position. Moreover, the noise level is reduced to a minimum.

A USB A-B cable for service purposes as well as connecting to ECG and pulmonary testing devices will be standard delivered with the product. To connect older ECG and pulmonary testing devices with RS232 or other connectors you need a special interface cable that can be ordered separately.

Features



Compatible with ECG and pulmonary devices

The Lode ergometers have both analog and digital interfaces and can be controlled easily by all known stress ECG and pulmonary devices available in the world. This is one of the reasons why the Lode ergometers are very popular worldwide.



Extreme low start up load

The extreme low start-up load of 7 watts and the adjustability in small steps of 1 watt make **att** this ergometer perfectly suitable for many different applications. The standard control unit shows multiple ergometry parameters and you can determine your specific default setting and start-up menu.



Low noise

Due to accurate manufacturing and the careful choice of materials the product has an extremely low noise level.



Accurate over a long period of time

The Lode ergometers are supplied with an electro-magnetic braking mechanism of Lanooy (eddy current). The biggest advantage of this braking system compared to a friction braking system is the absolute accuracy and the accuracy over time. Moreover, friction braking systems have more wearing parts.



RS232 connectivity

RS232 ports enable connectivity to most ECG and ergospirometry devices as well as PC's.



Readout out of saddle height

The height of the saddle is stepless adjustable and can be read-out on the saddle shaft



Perfect ergonomic position

Improved ergonomic position according to the latest requirements.



Ultra-low step-through

The lowest possible step-through guarantees easy access to the ergometer for all test subjects: a must for people who are not so mobile!



Hidden connectors

The cables are connected to the ergometer under the ergometer, which means that the test subject or operator cannot bump onto the connector.



USB connectivity

USB to connect to PC or ECG or ergospirometry products facilitates easy connectivity.



Can be controlled by all known Cardio Pulmonary Exercise devices



Corival cpet can a.o be extended with the following options:

Control Unit with touch screen 7" for ergometer Multifunctionality



Partnumber: 945834

Control Unit with touch screen 7" ordered addionally Multifunctionality



Partnumber: P945834

Programmable Control Unit with 7" Touchscreen Programmable



Partnumber: 945835

Blood Pressure Measurement with ECG trigger for biycle with ECG trigger



Partnumber: 945828

Electric adjustable saddle height


Easy and accurate positioning



Partnumber: 960810

0-Watt start-up system

Lowest possible startup power



Partnumber: 960805

Heart rate for bicycle ergometers

Heart rate in beats per minute



Partnumber: 945821

Adjustable cranks

Optimal force application



Partnumber: 928804

Shortened saddle shaft

Increase flexibility for smaller people



Partnumber: 960806

Saddle for children

Versatile ergometry



Partnumber: 401068

Saddle extra large

Versatile ergometry



Partnumber: 401084

Transportwheel for Corival

Easy transportation indoors



Partnumber: 960801

Network Module

Easy networking with LEM and LCRM



Partnumber: 945851

RS232 cable

Easy connection



Partnumber: 930911

USB to Serial converter

Easy connection



Partnumber: 226012



Can be controlled by all known Cardio Pulmonary Exercise devices



Specifications

Minimum load 7 W English user interface		
	•	
Maximum peak load 1000 W Norwegian user interface	~	
Minimum load increments 1 W Czech user interface	~	
Maximum continuous load 750 W Danish user interface	~	
Hyperbolic workload control ✓ Dutch user interface	~	
Linear workload control Finnish user interface	~	
Fixed torque workload control	~	
Maximum rpm independent constant load 150 rpm German user interface	~	
Minimum rpm independent constant load 30 rpm Italian user interface	~	
Optional heart rate controlled workload Japanese user interface	~	
Electromagnetic "eddy current" braking system	~	
Dynamic calibration Polish user interface	~	
Power range at maximum rpm (maximum) 1000 W Portugese user interface	~	
Accuracy Russian user interface	~	
Workload accuracy below 100 W 3 W Spanish user interface	~	
Workload accuracy from 100 to 500 W 3 % Turkish user interface	~	
Workload accuracy from 500 to 1000 W 5 % Ukrainian user interface	~	
Comfort Readout RPM	~	
Q-factor 180 mm Readout Time	~	
Minimum leg length user (incl. adjustable 602 mm 23.7 inch Readout Power pedals)	~	
Allowed user weight 180 kg 396.8 lbs Set Resistance	~	
Handlebar adjustment angle 360 ° Terminal operation mode	~	
Adjustability range seat 300 mm 11.8 inch Screen size (diagonal)	8.9 cm	3.5 inch
Minimum leg length user 645 mm 25.4 inch Touchscreen	~	
Connectivity		
Lode 38K4 interface protocol	~	
Lode interface protocol	~	
Lode WLP interface protocol	~	
Ergoline P10 interface protocol	~	
Ergoline P4 interface protocol	~	
Schiller interface protocol	~	
Bosch EKG 506 DS interface protocol	~	
USB connector	~	
RS232 in connector	~	
Dimensions		
Product length (cm)	105 cm	41.3 inch
Product width (cm)	46 cm	18.1 inch
Product height	114 cm	44.9 inch
Product weight	65 kg	143.3 lbs



Can be controlled by all known Cardio Pulmonary Exercise devices



Power requirements

Power cord length	250 cm	98.4 inch
100 - 240 V 50/60 Hz (160 Watt)	~	
Standards & Safety		
ISO 13485:2003 compliant	~	
ISO 9001:2008 compliant	~	
IEC 60601-1:2005	~	
Certification		
CE class Im according to MDD93/42/EEC	~	
CTüVus according to NRTL	~	
CB according to IECEE CB	~	

Order info

Partnumber:	960900	

^{*}Specifications are subject to change without notice.

